



VENEPUNCTURE ARMS

- *Standard*
- *Bag & Stand Advanced*
- *Advanced*

00330 00331

00440 00441

00290 00298



Learning has never been more lifelike

limbsandthings.com

STANDARD VENEPUNCTURE ARM

00330

00331

Simulated
Patient



Standard Venepuncture
Arm veins

BAG & STAND ADVANCED VENEPUNCTURE ARM

00440

00441

Simulated
Patient



Advanced Venepuncture
Arm veins



The Standard Venepuncture Arm provides the ideal platform for practising the skills related to venepuncture and IV cannulation.

The Bag & Stand Advanced Venepuncture Arm contains a more complex vein structure in combination with the simplicity of the bag and stand mock blood supply.

The bag and stand mock blood supply is perfect for a clinical skills training environment where large numbers of trainees are practising on a regular basis, or for self-directed learning time in a skills laboratory. Little technical knowledge is required for setup and maintenance.

VENEPUNCTURE ARMS

Product Qualities

Anatomy

- Soft, flexible skin
- Underlying palpable veins:
 - All models: cephalic, basilic and dorsal metacarpal
 - Advanced model only: median cubital

Key Features

- Realistic blood flashback
- Arm allows for digital pressure to stem blood flow
- Infusion tube allows for volume IV fluids to be administered

Cleaning & Care

- Veins are easy to replace
- Easy to clean, service and maintain

Skills Gained

- Venepuncture
- IV Cannulation
- IV Infusion
- Professional-to-patient communication

Versatility

- Skin is durable
- Veins last for up to 500 insertions of 21g needle
- Blood has been reformulated to improve leak resistance of veins (concentrated liquid to save you money and space)
- Can be used with vacuum blood collection systems, needle and syringe and IV cannulas
- Can be used for hybrid simulation and for professional-to-patient communication when using the optional Venepuncture Arm Harness

Safety

- Veins contain latex

ADVANCED VENEPUNCTURE ARM

00290

00298

Simulated Patient



Advanced Venepuncture Arm veins

NEW
SKIN
COLOUR

The Advanced Venepuncture Arm features a complex vein system, which is connected to a closed blood supply, for easy adjustment of pressure.

Contents

	STANDARD	BAG & STAND ADVANCED	ADVANCED
• 00292 Advanced Venepuncture Arm Veins	✗	✓	✓
• 00295 Advanced Venepuncture Vein Module	✗	✓	✓
• 00333 Standard Venepuncture Arm Veins	✓	✗	✗
• 00334 Standard Venepuncture Vein Module	✓	✗	✗
• 00335 Venepuncture Arm Infusion Tube	✓	✓	✓
• 00336 Bag & Stand Venepuncture Arm Mock Blood Supply	✓	✓	✗
• 00337 Pressurised Venepuncture Arm Mock Blood Supply	✗	✗	✓
• 00020 Concentrated Venous Blood - Starter Pack	✓	✓	✓
• 00302 Venepuncture Arm Rest	✓	✓	✓
• 00303 Refill Bottle	✗	✗	✓
• Venepuncture arm skin	✓	✓	✓
• Water-based lubricant pack (x3)	✓	✓	✓
• Arm shell	✓	✓	✓



Optional Extras & Consumables

(Available for all Venepuncture Arms)

- 00291 Arm Skin - Light
- 00299 Arm Skin - Brown
- 00300 Venepuncture Arm Harness
- 00301 Venepuncture Arm Carry Case



Curriculum References

Undergraduate Medicine

Appendix 1: Practical Procedure no. 4, p.12 – Venepuncture: Inserting a needle into a patient's vein to take a sample of blood for testing, or give an injection into the vein.

Source: *GMC Outcomes for Graduates: Tomorrow's Doctors*

Postgraduate Nursing

Follow guidelines for, undertake and record the following tasks: Venepuncture.

Source: *RCGP General Practice Foundation | General Practice Nurse competencies 2012*



For demonstration or to order

Limbs & Things Ltd. Sussex Street, St Phillips, Bristol, BS2 0RA, UK

T: 00 44 (0)117 311 0500 E: gbsales@limbsandthings.com

limbsandthings.com